

इंटरनेट

मानक

Disclosure to Promote the Right To Information

Whereas the Parliament of India has set out to provide a practical regime of right to information for citizens to secure access to information under the control of public authorities, in order to promote transparency and accountability in the working of every public authority, and whereas the attached publication of the Bureau of Indian Standards is of particular interest to the public, particularly disadvantaged communities and those engaged in the pursuit of education and knowledge, the attached public safety standard is made available to promote the timely dissemination of this information in an accurate manner to the public.

“जानने का अधिकार, जीने का अधिकार”

Mazdoor Kisan Shakti Sangathan

“The Right to Information, The Right to Live”

“पुराने को छोड़ नये के तरफ”

Jawaharlal Nehru

“Step Out From the Old to the New”

IS 4511-5 (2005): Methods of Test for Styrene-butadiene Rubber (SBR) Latices, Part 5: Determination of Soap and organic acid Content [PCD 13: Rubber and Rubber Products]



“ज्ञान से एक नये भारत का निर्माण”

Satyanarayan Gangaram Pitroda

“Invent a New India Using Knowledge”



“ज्ञान एक ऐसा खजाना है जो कभी चुराया नहीं जा सकता है”

Bhartrhari—Nitiśatakam

“Knowledge is such a treasure which cannot be stolen”

BLANK PAGE



भारतीय मानक

स्टाइरीन-ब्यूटाडाइन रबड़ (एसबीआर),

लेटेसिस की परीक्षण पद्धति

भाग 5 साबुन एवं कार्बनिक-अम्ल का अंश ज्ञात करना

(दूसरा पुनरीक्षण)

Indian Standard

**METHOD OF TEST FOR STYRENE-BUTADIENE
RUBBER (SBR), LATICES**

PART 5 DETERMINATION OF SOAP AND ORGANIC-ACID CONTENT

(Second Revision)

ICS 83.140.40

© BIS 2005

BUREAU OF INDIAN STANDARDS
MANAK BHAVAN, 9 BAHADUR SHAH ZAFAR MARG
NEW DELHI 110002

NATIONAL FOREWORD

This Indian Standard (Part 5) (Second Revision) which is identical with ISO 7781 : 2001 'Styrene-butadiene rubber, raw — Determination of soap and organic-acid content' issued by the International Organization for Standardization (ISO) was adopted by the Bureau of Indian Standards on the recommendations of the Rubber and Rubber Products Sectional Committee and approval of the Petroleum, Coal and Related Products Division Council.

The text of ISO Standard has been proposed to be approved as suitable for publication as an Indian Standard without deviations. Certain conventions are, however, not identical to those used in Indian Standards. Attention is particularly drawn to the following:

- a) Wherever the words 'International Standard' appear referring to this standard, they should be read as 'Indian Standard'.
- b) Comma (,) has been used as a decimal marker while in Indian Standards, the current practice is to use a point (.) as the decimal marker.

The Technical Committee responsible for the preparation of this standard has reviewed the provisions of the following International Standards and has decided that they are acceptable for use in conjunction with this standard:

<i>International Standard</i>	<i>Title</i>
ISO 385-1 : 1984	Laboratory glassware — Burettes — Part 1 : General requirements
ISO 648 : 1977	Laboratory glassware — One-mark pipettes
ISO 1042 : 1998	Laboratory glassware — One-mark volumetric flasks
ISO 1795 : 2000	Rubber, raw natural and raw synthetic — Sampling and further preparative procedures
ISO 4799 : 1978	Laboratory glassware — Condensers

In reporting the results of a test or analysis made in accordance with this standard, if the final value, observed or calculated, is to be rounded off, it shall be done in accordance with IS 2 : 1960 'Rules for rounding off numerical values (*revised*)'.

Indian Standard

METHOD OF TEST FOR STYRENE-BUTADIENE RUBBER (SBR), LATICES

PART 5 DETERMINATION OF SOAP AND ORGANIC-ACID CONTENT

(*Second Revision*)

WARNING — Persons using this International Standard should be familiar with normal laboratory practice. This standard does not purport to address all of the safety problems, if any, associated with its use. It is the responsibility of the user to establish appropriate safety and health practices and to ensure compliance with any national regulatory conditions.

1 Scope

This International Standard specifies methods for the determination of the soap and organic-acid content of raw styrene-butadiene rubber (SBR). Method A is suitable for non-oil-extended styrene-butadiene rubbers. Method B is applicable to oil-extended styrene-butadiene rubbers.

The methods depend on the extraction of the organic acids and soaps from the rubber by means of a specified solvent. In practice, therefore, it is convenient to determine both organic-acid and soap contents on separate portions of the same solvent extract. Since the soaps and organic acids present in the rubber are not single chemical compounds, the method gives only an approximate value for the soap and organic-acid content.

2 Normative references

The following normative documents contain provisions which, through reference in this text, constitute provisions of this International Standard. For dated references, subsequent amendments to, or revisions of, any of these publications do not apply. However, parties to agreements based on this International Standard are encouraged to investigate the possibility of applying the most recent editions of the normative documents indicated below. For undated references, the latest edition of the normative document referred to applies. Members of ISO and IEC maintain registers of currently valid International Standards.

ISO 385-1:1984, *Laboratory glassware — Burettes — Part 1: General requirements*

ISO 648:1977, *Laboratory glassware — One-mark pipettes*

ISO 1042:1998, *Laboratory glassware — One-mark volumetric flasks*

ISO 1795:2000, *Rubber, raw natural and raw synthetic — Sampling and further preparative procedures*

ISO 4799:1978, *Laboratory glassware — Condensers*

3 Principle

A weighed test portion of the rubber, in the form of thin strips, is extracted by ethanol-toluene azeotrope, or, for alum-coagulated rubber, ethanol-toluene-water mixture. After making up to a standard volume, an aliquot portion of the extract is withdrawn and titrated with standard acid for the determination of soap and with standard alkali for the determination of organic acid.

4 Reagents

During the analysis, use only reagents of recognized analytical grade and only distilled water or water of equivalent purity. Use commercially available analytical-grade buffer solutions of known pH or, in the absence of commercial buffer solutions, prepare the solutions required (4.6, 4.7 and 4.8).

4.1 Ethanol-toluene azeotrope (ETA).

Mix 7 volumes of absolute ethanol with 3 volumes of toluene. Alternatively, mix 7 volumes of commercial-grade ethanol with 3 volumes of toluene, and boil the mixture with anhydrous calcium oxide (quicklime) under reflux for 4 h. Cool to room temperature and decant through No. 42 filter paper.

4.2 Ethanol-toluene-water mixture.

Mix 95 cm³ of ETA (4.1) and 5 cm³ of water.

4.3 Sodium hydroxide solution, $c(\text{NaOH}) = 0,1 \text{ mol/dm}^3$, accurately standardized.

4.4 Thymol blue indicator.

Dissolve 0,06 g of thymol blue in 6,45 cm³ of 0,02 mol/dm³ sodium hydroxide solution and dilute to 50 cm³ with water.

4.5 Hydrochloric acid, $c(\text{HCl}) = 0,05 \text{ mol/dm}^3$, accurately standardized.

4.6 Buffer solution of nominal pH 7.

Dissolve 3,40 g of potassium dihydrogen phosphate (KH_2PO_4) and 3,55 g of disodium hydrogen phosphate (Na_2HPO_4) in water and make up to 1 000 cm³ in a volumetric flask. The pH of this solution is 6,87 at 23 °C.

Store the solution in a glass or polyethylene vessel that is resistant to chemicals.

4.7 Buffer solution of nominal pH 4.

Dissolve 10,21 g of potassium hydrogen phthalate ($\text{HOOC}_6\text{H}_4\text{COOK}$) in water and make up to 1 000 cm³ in a volumetric flask. The pH of this solution is 4,00 at 23 °C.

Store the solution in a glass or polyethylene vessel that is resistant to chemicals.

4.8 Buffer solution of nominal pH 9.

Dissolve 3,814 g of sodium tetraborate decahydrate ($\text{Na}_2\text{B}_4\text{O}_7 \cdot 10\text{H}_2\text{O}$) in water and make up to 1 000 cm³ in a volumetric flask. The pH of this solution, when freshly prepared, is 9,20 at 23 °C.

Store the solution in a glass or polyethylene vessel that is resistant to chemicals and fitted with a soda lime carbon dioxide trap. Replace the solution after 1 month.

5 Apparatus

5.1 Balance, accurate to 1 mg.

5.2 Hotplate.

5.3 Wide-mouthed conical flask, capacity 400 cm³ to 500 cm³.

5.4 Volumetric flask, capacity 250 cm³, complying with the requirements of ISO 1042.

5.5 Reflux condenser, complying with the requirements of ISO 4799.

5.6 Conical flask, capacity 250 cm³.

NOTE Alternatively, a Soxhlet extractor may be used instead of a reflux condenser and a conical flask.

5.7 Burette, capacity 25 cm³, complying with the requirements of ISO 385-1.

5.8 Pipette, capacity 100 cm³, complying with the requirements of ISO 648.

5.9 Automatic titrator or pH-meter, including a calomel electrode and glass electrode, accurate to within 10 mV or 0,1 pH-units.

5.10 Magnetic stirrer, with a polytetrafluoroethylene-coated stirrer bar.

6 Sampling and sample preparation

Sheet out 2 g to 6 g of rubber, selected and prepared in accordance with ISO 1795. Cut into pieces no larger than 2 mm × 2 mm or strips no longer than 10 mm or wider than 5 mm. Weigh a test portion of approximately 2 g to the nearest 0,001 g.

7 Procedure

7.1 Preparation of the test solution

Place a circular filter paper in the bottom of the wide-mouthed conical flask (5.3) and add 100 cm³ of ETA extraction solvent (4.1) for all rubbers except alum-coagulated rubbers. For alum-coagulated rubbers, use ethanol-toluene-water mixture (4.2).

Introduce the strips of rubber separately into the flask, swirling after each addition so that the strips are thoroughly wetted with solvent and sticking is minimized.

Fit the reflux condenser (5.5) to the flask and boil the solvent very gently under reflux for 1 h.

Transfer the extract to the volumetric flask (5.4), and treat the rubber with a second 100 cm³ portion of extraction solvent under reflux for 1 h. Add this extract also to the volumetric flask. Rinse the strips with three successive 10 cm³ portions of extraction solvent, add these washings to the volumetric flask and, after cooling to room temperature, adjust the final volume to 250 cm³ with appropriate extraction solvent.

NOTE Alternatively, the weighed strips of sample may be wrapped in filter paper and placed in a Soxhlet extractor (see note to 5.6) and extracted with ETA (4.1) or ethanol-toluene-water mixture (4.2) under reflux for a minimum of 4 h.

7.2 Method A (for non-oil-extended SBR)

7.2.1 Procedure for determination of soap content

After thorough mixing, pipette 100 cm³ of the diluted extract into the 250 cm³ conical flask (5.6), add six drops of thymol blue indicator (4.4) and titrate the solution with hydrochloric acid solution (4.5) to the first colour change (V_1).

Carry out a blank titration on 100 cm³ of extraction solvent taken from the same stock as was used for the test portion and using the same indicator as was used for titration of the test portion (V_2).

7.2.2 Procedure for determination of organic-acid content

Proceed exactly as in 7.2.1, but titrate the aliquot portion with sodium hydroxide solution (4.3), again using thymol blue indicator (4.4) (V_3).

Carry out a blank titration on 100 cm³ of extraction solvent using the same method (V_4).

7.3 Method B (for oil-extended SBR)

7.3.1 Procedure for determination of soap content

Turn on the autotitrator or pH-meter (5.9), and allow the electronic circuit to stabilize. The temperature shown by the temperature compensator indicator shall be same as the temperature of the test solution.

Calibrate the autotitrator or pH-meter using nominally pH 7 buffer solution (4.6) and pH 4 buffer solution (4.7).

Pipette 100 cm³ of the test solution into a 250 cm³ beaker containing a stirrer bar, then place the beaker on the magnetic stirrer (5.10). Insert a glass electrode and a calomel electrode into the test solution in the beaker. While stirring, titrate the test solution with hydrochloric acid solution (4.5) to pH 4,8, slowing the titration rate near the equivalence point. Record the volume of hydrochloric acid solution used at the equivalence point (V_1).

Carry out a blank titration on 100 cm³ of extraction solvent using the same method (V_2).

7.3.2 Procedure for determination of organic-acid content

Calibrate the autotitrator or pH-meter using nominally pH 7 buffer solution (4.6) and pH 9 buffer solution (4.8).

Proceed exactly as in 7.3.1, but titrate the aliquot portion with sodium hydroxide solution (4.3) to pH 11,5, again slowing the titration rate near the equivalence point. Record the volume of sodium hydroxide solution used at the equivalence point (V_3).

Carry out a blank titration on 100 cm³ of extraction solvent using the same method (V_4).

8 Expression of results

8.1 Calculate the soap content using the equation

$$w_s = \frac{0,25 \times (V_1 - V_2) \times c_1 \times K_s}{m}$$

where

w_s is the soap content, as a percentage by mass;

V_1 is the volume, in cubic centimetres, of hydrochloric acid solution used to titrate the rubber extract;

V_2 is the volume, in cubic centimetres, of hydrochloric acid solution used to titrate the blank;

c_1 is the actual concentration, in mol/dm³, of the hydrochloric acid solution (4.5);

m is the mass, in grams, of the test portion;

K_s is the appropriate factor selected from the following:

306 when the soap is to be calculated as sodium stearate,

368 when the soap is to be calculated as sodium rosinate,

337 when the soap is to be calculated as a 50:50 mixture of sodium stearate and sodium rosinate,

322 when the soap is to be calculated as potassium stearate,

384 when the soap is to be calculated as potassium rosinate,

353 when the soap is to be calculated as a 50:50 mixture of potassium stearate and potassium rosinate,

345 when the soap is to be calculated as a 50:50 mixture of sodium stearate and potassium rosinate or of sodium rosinate and potassium stearate.

NOTE Since the soaps present in the rubber are not single chemical compounds, the value assigned to K_s gives only an approximate value for the soap content. A test for rosin is given in annex A.

8.2 Calculate the organic-acid content using the equation

$$w_0 = \frac{0,25 \times (V_3 - V_4) \times c_2 \times K_0}{m}$$

where

w_0 is the organic-acid content, expressed as a percentage by mass;

V_3 is the volume, in cubic centimetres, of sodium hydroxide solution used to titrate the test solution;

V_4 is the volume, in cubic centimetres, of sodium hydroxide solution used to titrate the blank;

c_2 is the actual concentration, in mol/dm³, of the sodium hydroxide solution (4.3);

m is the mass, in grams, of the test portion;

K_0 is the appropriate factor selected from the following:

284 when the acid is to be calculated as stearic acid,

346 when the acid is to be calculated as rosin acid,

315 when the acid is to be calculated as a 50:50 mixture of stearic acid and rosin acid.

NOTE Since the organic acids present in the rubber are not single chemical compounds, the value assigned to K_0 gives only an approximate value for the organic acid content. A test for rosin is given in annex A.

9 Test report

The test report shall include the following information:

- all details required for the identification of the sample;
- a reference to this International Standard;
- the soap and/or organic-acid content and the value of K_0 used;
- the method used (A or B);
- details of any deviation from this International Standard;
- the date of the test.

Annex A **(informative)**

Test for rosin

A.1 Reagents

A.1.1 Acetic anhydride.

A.1.2 Sulfuric acid solution.

Carefully add 65 g of sulfuric acid ($\rho_{20} = 1,84 \text{ g/cm}^3$) to 35 g of water and mix well.

A.1.3 Potassium permanganate solution, $c(\text{KMnO}_4) = 0,0002 \text{ mol/dm}^3$.

A.2 Procedure

Mix a small amount of sample with about 3 cm^3 of acetic anhydride (A.1.1). Add 2 drops of sulfuric acid (A.1.2). The reaction is positive for rosin if a temporary violet colour appears which, at the moment of its maximum intensity, is stronger than the colour of the potassium permanganate solution.

Bureau of Indian Standards

BIS is a statutory institution established under the *Bureau of Indian Standards Act, 1986* to promote harmonious development of the activities of standardization, marking and quality certification of goods and attending to connected matters in the country.

Copyright

BIS has the copyright of all its publications. No part of these publications may be reproduced in any form without the prior permission in writing of BIS. This does not preclude the free use, in the course of implementing the standard, of necessary details, such as symbols and sizes, type or grade designations. Enquiries relating to copyright be addressed to the Director (Publications), BIS.

Review of Indian Standards

Amendments are issued to standards as the need arises on the basis of comments. Standards are also reviewed periodically; a standard along with amendments is reaffirmed when such review indicates that no changes are needed; if the review indicates that changes are needed, it is taken up for revision. Users of Indian Standards should ascertain that they are in possession of the latest amendments or edition by referring to the latest issue of 'BIS Catalogue' and 'Standards: Monthly Additions'.

This Indian Standard has been developed from Doc: No. PCD 13 (2168).

Amendments Issued Since Publication

Amend No.	Date of Issue	Text Affected

BUREAU OF INDIAN STANDARDS

Headquarters:

Manak Bhavan, 9 Bahadur Shah Zafar Marg, New Delhi 110002

Telephones: 2323 0131, 2323 3375, 2323 9402

website : www.bis.org.in

Regional Offices:

	Telephones
Central : Manak Bhavan, 9 Bahadur Shah Zafar Marg NEW DELHI 110002	{ 2323 7617 2323 3841
Eastern : 1/14 C.I.T. Scheme VII M, V.I.P. Road, Kankurgachi KOLKATA 700054	{ 2337 8499, 2337 8561 2337 8626, 2337 9120
Northern : SCO 335-336, Sector 34-A, CHANDIGARH 160022	{ 260 3843 260 9285
Southern : C.I.T. Campus, IV Cross Road, CHENNAI 600113	{ 2254 1216, 2254 1442 2254 2519, 2254 2315
Western : Manakalaya, E9 MIDC, Marol, Andheri (East) MUMBAI 400093	{ 2832 9295, 2832 7858 2832 7891, 2832 7892

Branches: AHMEDABAD. BANGALORE. BHOPAL. BHUBANESHWAR. COIMBATORE. FARIDABAD. GHAZIABAD. GUWAHATI. HYDERABAD. JAIPUR. KANPUR. LUCKNOW. NAGPUR. NALAGARH. PATNA. PUNE. RAJKOT. THIRUVANANTHAPURAM. VISAKHAPATNAM.